



Anti-Mouse IgG polyclonal antibody [HRP] (CPBT-68096GM)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat anti Mouse IgG (H/L)Goat anti Mouse IgG (H/L) antibody recognises mouse IgG heavy and light chains.
Specificity	IgG
Immunogen	Native murine immunoglobulin G.
Isotype	IgG
Source/Host	Goat
Species Reactivity	Mouse
Conjugate	HRP
Applications	IHC-Fr; ELISA; WB
Format	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Size	2 ml
Preservative	None
Storage	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

Gene Name	Ighg1 immunoglobulin heavy constant gamma 1 (G1m marker) [Mus musculus (house mouse)]
------------------	---

Official Symbol	IGHG1
Synonyms	IGHG1; immunoglobulin heavy constant gamma 1 (G1m marker); IgG1; Igh-4; VH7183; IGG;
Entrez Gene ID	16017
UniProt ID	P01837
Chromosome Location	12 F1; 12 62.1 cM
Function	antigen binding;