



Anti-Mouse IgG polyclonal antibody [AP] (CPBT-68110GM)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	This product recognises Mouse IgG. Cross-reactivity with Bovine, Chicken, Goat, Guinea pig, Hamster, Horse, Human, Rabbit, Rat or Sheep immunoglobulin has been minimised by solid phase absorption.
Specificity	IgG
Immunogen	Mouse IgG whole molecule
Isotype	IgG
Source/Host	Goat
Species Reactivity	Mouse
Conjugate	AP
Applications	IHC-Fr; ELISA; WB
Format	Purified IgG conjugated to Alkaline Phosphatase - liquid
Size	1 mg
Preservative	0.1% Sodium Azide
Storage	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted.

GENE INFORMATION

Gene Name	Ighg1 immunoglobulin heavy constant gamma 1 (G1m marker) [Mus musculus (house mouse)]
------------------	---

Official Symbol	IGHG1
Synonyms	IGHG1; immunoglobulin heavy constant gamma 1 (G1m marker); IgG1; Igh-4; VH7183; IGG;
Entrez Gene ID	16017
UniProt ID	P01837
Chromosome Location	12 F1; 12 62.1 cM
Function	antigen binding;