



Anti-Mouse IgG + IgM + IgA polyclonal antibody [FITC] (CPBT-68125GM)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat anti Mouse IgG/A/M (Human Adsorbed) antibody recognizes all subclasses of murine immunoglobulin and has been adsorbed against human, bovine and hamster immunoglobulins. Flow Cytometry Use 10ul of the suggested working dilution to label 106 cells in 100ul.
Specificity	IgG IgA IgM
Isotype	IgG
Source/Host	Goat
Species Reactivity	Mouse
Conjugate	FITC
Applications	FC
Format	Purified IgG - liquid
Preservative	0.09% Sodium Azide
Storage	in frost free freezers is not recommended. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

BACKGROUND

Introduction	Immunoglobulin (Ig) is a large Y-shape protein produced by plasma cells that is used by the immune system to identify and neutralize foreign objects such as bacteria and viruses. The antibody recognizes a unique part of the foreign target, called an antigen.
---------------------	--

Keywords

IgG; IgA; IgM;
