



Anti-HISTONE H3 ACETYLATED polyclonal antibody (CPBT-68206RH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Histone H3 acetylated at lysine 14 (will recognise all or most acetylated H3 isoforms).
Specificity	HISTONE H3 ACETYLATED
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Amphibia, Mammals, Yeast
Conjugate	Unconjugated
Applications	ELISA; IF; WB
Format	Serum - liquid
Size	100 µl
Preservative	0.09% Sodium Azide
Storage	in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

Gene Name	H3F3C H3 histone, family 3C [Homo sapiens (human)]
Official Symbol	H3F3C

Synonyms	H3F3C; H3 histone, family 3C; H3.5; histone H3.3C; histone H3.5; histone variant H3.5; HISTONE H3 ACETYLATED;
Entrez Gene ID	440093
Protein Refseq	NP_001013721
UniProt ID	Q6NXT2
Chromosome Location	12p11.21
Pathway	Alcoholism; Systemic lupus erythematosus; Transcriptional misregulation in cancer;
Function	nucleosomal DNA binding; protein heterodimerization activity;