



Recombinant HIV P17 [His] (DAG-2799)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------|---|
| Antigen Description | HIV p17 is the first HIV antigen, which is developed for some high quality HIV diagnostic kits. |
| Species | HIV |
| Purity | >95% pure, SDS Gel Electrophoresis Result |
| Conjugate | His |
| Molecular Weight | p17~29kD |
| Format | Purified, Liquid |
| Concentration | 5.0 mg/ml |
| Size | 200 µg |
| Buffer | PBS, pH 7.4±0.1 |
| Preservative | None |
| Storage | Store at 4 °C. Avoid multiple freeze/thaw cycles. |

BACKGROUND

| | |
|----------|--|
| Keywords | HIV p17; human immunodeficiency virus p17; human immunodeficiency virus; HIV |
|----------|--|