



HIV-M2 (gp41+gp120+gp36+O) (DAG-B910)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Expression System	E. coli
Purity	Purity >90%
Conjugate	Unconjugated
Procedure	None
Format	Purified, protein solution in Tris-HCl pH 10.0
Size	1 mg
Buffer	70mM Tris-HCl pH 10.0
Preservative	None
Storage	Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles.
Ship	IA