



## E. coli BirA [His] (DAG-H10051)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|              |  |
|--------------|--|
| Species      | E. coli  |
| Purity       | > 95 % as determined by SDS-PAGE   |
| Conjugate    | His  |
| Size         | 50 µg, 100 µg  |
| Preservative | None   |
| Storage      | Store it under sterile conditions at -70 °C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. |

### BACKGROUND

|              |  |
|--------------|--|
| Introduction | No Information.  |
| Keywords     | birA; biotin ligase; Bifunctional protein birA; Bio R; BioR; Biotin [acetyl CoA carboxylase] synthetase; Biotin operon repressor; Biotin protein ligase; Bir A; BirA bifunctional protein; Dhb B; dhbB; biotin--protein ligase |