



HSV Grade 2 (DAG-H10363)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Grade 2 antigen contains a high concentration of virus particles and components in tissue culture medium containing a low concentration of bovine serum. The preparation contains no preservative or detergent. Some residual host cellular material is present
Species	HSV
Conjugate	Unconjugated
Applications	gG and IgM ELISA Western Blot Lymphocyte proliferation assays
Size	1 ml, 10 ml
Preservative	None
Storage	Store this antigen preparation frozen at -70 °C to -1 00 °C. Repeated freezing and thawing should be avoided

BACKGROUND

Introduction	Herpes simplex virus 1 and 2 (HSV-1 and HSV-2), also known as Human herpes virus 1 and 2 (HHV-1 and -2), are two members of the herpes virus family, Herpesviridae, that infect humans.[1] Both HSV-1 (which produces most cold sores) and HSV-2 (which produce
Keywords	Herpes simplex virus 2; HSV-2; Human herpes virus 2; HSV2