



E. coli Gor (aa 1 - 450) (DAG-P2176)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	E. coli gor full length protein
Species	E. coli
Purity	> 90 % by SDS-PAGE. This antigen is purified using conventional chromatography techniques.
Conjugate	Unconjugated
Applications	SDS-PAGE
Sequence Similarities	Belongs to the class-I pyridine nucleotide-disulfide oxidoreductase family.
Molecular Weight	51 kDa including tags
Format	Liquid
Buffer	pH: 8.00 Constituents: 0.02% DTT, 0.32% Tris HCl, 10% Glycerol, 0.58% Sodium chloride
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. pH: 8.00 Constituents: 0.02% DTT, 0.32% Tris HCl, 10% Glycerol, 0.58% Sodium chloride

BACKGROUND

Introduction	Escherichia coli; commonly abbreviated E. coli) is a gram-negative, facultatively anaerobic, rod-shaped bacterium of the genus Escherichia that is commonly found in the lower intestine of warm-blooded organisms (endotherms). Most E. coli strains are harml
Keywords	ECK3485; glutathione oxidoreductase; Glutathione reductase; gor; gorA; GR; Grase; GSHR_ECOLI; JW3467; E. coli gor; Escherichia coli gor