



HSV type 2 (DAG-P2201)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	HSV HSV2 full length protein
Species	HSV
Conjugate	Unconjugated
Format	Liquid
Buffer	Constituents: 99% PBS, 0.1% Triton-X-100
Preservative	None
Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. Constituents: 99% PBS, 0.1% Triton-X-100

BACKGROUND

Introduction	Herpes simplex virus 1 and 2 (HSV-1 and HSV-2), also known as human herpesvirus 1 and 2 (HHV-1 and HHV-2), are two members of the herpesvirus family, Herpesviridae, that infect humans.
Keywords	Herpes Simplex Virus Type 2; HSV 2