



# HIV type 1 P311 (full length) (DAG-P2466)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	HIV-1 p31 full length protein
<b>Species</b>	HIV
<b>Purity</b>	>95% by SDS-PAGE .
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA WB
<b>Format</b>	Liquid
<b>Buffer</b>	Information available upon request.
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Store at -20°C. Information available upon request.

## BACKGROUND

<b>Introduction</b>	The human immunodeficiency virus (HIV) is a lentivirus (slowly replicating retrovirus) that causes the acquired immunodeficiency syndrome (AIDS), a condition in humans in which progressive failure of the immune system allows life-threatening opportunistic
<b>Keywords</b>	HIV1 p311