



Carbofuran [BSA] (DAG-P2540B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Carbofuran is one of the most toxic carbamate pesticides, also named Furdan in the market. It is used to control insects in a wide variety of field crops, including potatoes, corn and soybeans. Since it is a systemic insecticide, the plant absorbs it through the roots, and from here the plant distributes it throughout its organs where insecticidal concentrations are attained. Carbofuran also has contact activity against pests. But Carbofuran is banned in Canada and the European Union nowadays.
Species	N/A
Conjugate	BSA
Applications	The antigen is intended for use in ELISA assay.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.