



Recombinant HIV type 1 P24 Protein [GST] (DAG-P2664)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	HIV1 p24 full length protein
Antigen Description	Capsid protein p24 probably forms the conical core of the virus that encapsulates the genomic RNA-nucleocapsid complex.
Species	HIV
Purity	> 90 % by SDS-PAGE.
Conjugate	GST
Applications	ELISA WB
Bio-activity	Strongly reacts with human HIV positive serum.
Format	Liquid
Buffer	Preservative: None Constituents: 50% Glycerol, 1.5M Urea, 25mM Tris HCl, pH 8.0
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. Preservative: None Constituents: 50% Glycerol, 1.5M Urea, 25mM Tris HCl, pH 8.0

BACKGROUND

Introduction	The human immunodeficiency virus (HIV) is a lentivirus (slowly replicating retrovirus) that causes the acquired immunodeficiency syndrome (AIDS), a condition in humans in which progressive failure of the immune system allows life-threatening opportunistic
---------------------	---

Keywords

CA; Capsid protein p24; Human immunodeficiency virus 1; Human immunodeficiency virus type 1 p24
