



Recombinant EBV EBNA1 mosaic Protein (a.a. 1-90, 408-498) (DAG-T1130)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	E.coli derived recombinant. The mosaic protein contains the HHV-4 EBNA regions (1-90/408-498aa).
Species	EBV
Conjugate	Unconjugated
Preservative	None
Storage	Store at 2–8 °C

BACKGROUND

Introduction	Epstein Barr virus (EBV) is a member of the herpes virus family and one of the most common human viruses. Most people become infected with EBV during their lives. Primary infections usually results in infectious mononucleosis (glandular fever) but the virus can also lay dormant in B lymphocytes and when reactivated become associated with more serious disease such as Burkitt's lymphoma, nasopharyngeal carcinoma and Hodgkin's disease.
Keywords	Epstain-Barr Virus EBNA1; EBV EBNA1; Epstain-Barr Virus; EBV