



GBM Antigen (DAG-T1242)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Purified from bovine kidney. After coating onto ELISA plates the product will bind autoantibodies to GBM antigen.
Purity	The GBM antigen (NC1 domain of type IV collagen from glomerular basement membrane) is more than 90% pure, as assessed by SDS polyacrylamide gel electrophoresis.
Procedure	None
Concentration	0.1-2.0 mg/ml
Size	0.20 mg
Preservative	None
Storage	The product is stabilised with 0.1% Micr-O-protect TM. Store at -20°C or below (long term) or at +4°C (short term). Avoid repeated freezing and thawing. Mix thoroughly before use.

BACKGROUND

Keywords	GBM;GBM Antigen;
-----------------	------------------