



Recombinant HIV Glycoprotein 41 [Gold] (DAG-T1426)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	40 nm particle coated with rHIV gp120/gp41 antigen at 10 micrograms per mL
Species	HIV
Conjugate	Gold
Size	1 ml, 10 ml, 100 ml
Buffer	Bovine Serum Albumin/Phosphate pH 7.4
Preservative	0.1% Sodium Azide
Storage	2°C - 8°C