



DENV type 1 protein (DAG-T2660)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Dengue fever is a viral disease that affects people in the tropics and certain regions of Africa. Related to malaria, it is caused by four different version of the flavivirus, and is also known as "breakbone fever."
Conjugate	Unconjugated
Applications	ELISA, IA, WB
Size	1 ml
Buffer	Supplied in PBS buffer
Preservative	None
Storage	Store at -70°C to -100°C. Avoid freeze / thaw cycles

BACKGROUND

Keywords	Dengue Type 1; Dengue virus Type 1 protein; Dengue virus protein; DENV protein
-----------------	--