



Aluminum-EDTA [BSA] (DAG-WT002)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Aluminum-EDTA conjugated with BSA
Specificity	This product is ELISA-paired with CABT-Z520M
Species	N/A
Conjugate	BSA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Aluminum powder, uncoated appears as a light gray solid. Denser than water. Contact may burn skin, eyes, and mucous membranes. May be toxic by ingestion, inhalation and skin absorption. Used to make other chemicals.
--------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------