



Tebuconazole [BSA] (DAG-WT022)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Tebuconazole conjugated with BSA
Specificity	This product is ELISA-paired with CABT-Z530M
Species	N/A
Conjugate	BSA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Tebuconazole is a triazole fungicide used agriculturally to treat plant pathogenic fungi. Its acute toxicity is moderate. Its bactericidal mechanism is mainly to inhibit the biosynthesis of ergosterol of pathogenic bacteria, diseases caused by Puccinia spp, Ceratocystis spp, Sclerotium spp and Needle spores.
--------------	---