



Penconazole [OVA] (DAG-WT023)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Penconazole conjugated with OVA
Specificity	This product is ELISA-paired with CABT-Z531M
Nature	Synthetic
Expression System	N/A
Species	N/A
Conjugate	OVA
Procedure	None
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Penconazole is a typical triazole chiral fungicide. It is widely used to control the pathogenic bacteria of powdery mildew, black genus, basidiomycetes, and deuteromycetes.
---------------------	--