



## Triazamidine [OVA] (DAG-WT025)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Product Overview</b>	Triazamidine conjugated with OVA
<b>Specificity</b>	This product is ELISA-paired with CABT-Z532M
<b>Species</b>	N/A
<b>Conjugate</b>	OVA
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M PBS, pH7.4
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

<b>Introduction</b>	Triazamidine belongs to the class of aromatic diamidines. It is a traditionally used broad-spectrum anti-blood protozoan. For trypanosomiasis caused by trypanosomiasis and equine disease. This product selectively blocks the DNA synthesis and replication of the trypanosomes' motifs, and irreversibly combines with the nucleus, so that the trypanosomes' motifs disappear and cannot divide and reproduce.
---------------------	--