



Pyridaben [BSA] (DAG-WT028)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Pyridaben conjugated with BSA
Specificity	This product is ELISA-paired with CABT-Z533M
Species	N/A
Conjugate	BSA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Pyridaben is a pyridazinone, an organochlorine insecticide and an organochlorine acaricide. It has a role as a mitochondrial NADH:ubiquinone reductase inhibitor.
--------------	---