



Butralin [BSA] (DAG-WT044)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Butralin conjugated with BSA
Specificity	This product is ELISA-paired with CABT-Z541M
Nature	Synthetic
Expression System	N/A
Species	N/A
Conjugate	BSA
Procedure	None
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Butraline is a racemate from the group of dinitroaniline derivatives, which is used as a pre-emergence herbicide. Butraline can be obtained by chlorinating tert-butylbenzene and then reacting it with nitric acid and sec-butylamine.
---------------------	---