



Pendimethalin [BSA] (DAG-WT046)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Pendimethalin conjugated with BSA
Specificity	This product is ELISA-paired with CABT-Z542M
Nature	Synthetic
Expression System	N/A
Species	N/A
Conjugate	BSA
Procedure	None
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Pendimethalin is a member of the class of substituted anilines that is N-(pentan-3-yl)aniline bearing two additional nitro substituents at positions 2 and 6 as well as two methyl substituents at positions 3 and 4. A herbicide used to control most annual grasses and many annual broad-leaved weeds. It has a role as a herbicide, an environmental contaminant and an agrochemical.
---------------------	---