



Triiodothyronine (T3) [OVA] (DAG-WT100)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Triiodothyronine, OVA-Conjugated
Conjugate	OVA
Applications	Immunoassay
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	The thyroid hormones, thyroxine (T4) and triiodothyronine (T3), are tyrosine-based hormones produced by the thyroid gland primarily responsible for regulation of metabolism. An important component in the synthesis of thyroid hormones is iodine. The major form of thyroid hormone in the blood is thyroxine (T4), which has a longer half life than T3.
--------------	--