



Triiodothyronine (T3) [AP] (DAG-WT230)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	AP
Applications	Immunoassays
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	50 Units
Buffer	0.06M Tris saline buffer with 0.1% BSA, 0.1% NaN ₃ and 0.1% MgCl
Preservative	None
Storage	Store at 4 °C.

BACKGROUND

Introduction	Triiodothyronine, C ₁₅ H ₁₂ I ₃ NO ₄ , also known as T ₃ , is a thyroid hormone. It affects almost every physiological process in the body, including growth and development, metabolism, body temperature, and heart rate.
Keywords	T ₃ ; 3,3',5-Triiodo-L-thyronine; Triiodothyronine; 3,5,3'-Triiodothyronine