



Recombinant Mus musculus Mus m 1 (DAG-WT279)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Nature	Recombinant
Expression System	E.coli
Purity	> 90%
Conjugate	TBD
Applications	Immunoassays
Molecular Weight	17 kDa
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	20mM Tris, 150mM NaCl, pH8.0
Preservative	None
Storage	Store at -80°C. Avoid freeze / thaw cycle
Warnings	PLEASE note that this product is intended for research use only; not for diagnostic or clinical use.

BACKGROUND

Keywords Mus musculus allergen 1; Mus m 1; Lipocalin and urinary prealbumin; Mus m 1.0101; Mus m

