



Recombinant Prunus persica Pru p 1 (DAG-WT288)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Nature	Recombinant
Expression System	E.coli
Purity	> 90%
Conjugate	TBD
Applications	Immunoassays
Molecular Weight	18.7 kDa
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	100mMPBS, pH7.4, 10% Glycerol
Preservative	None
Storage	Store at -80°C. Avoid freeze / thaw cycle
Warnings	PLEASE note that this product is intended for research use only; not for diagnostic or clinical use.

BACKGROUND

Keywords	Prunus persica allergen 1; Pru p 1.0101; Pru p 1.0201; Pru p 1.0301
-----------------	---