



Recombinant *Metapenaeus ensis* Met e 1 (DAG-WT293)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90%
Conjugate	TBD
Applications	Immunoassays
Molecular Weight	33 kDa
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	20mM Tris, 10mM NaCl, pH8.0
Preservative	None
Storage	Store at -80°C. Avoid freeze / thaw cycle

BACKGROUND

Keywords	<i>Metapenaeus ensis</i> allergen 1; Met e 1; Tropomyosin; Met e 1.0101
----------	---