



# Recombinant *Gadus callarias* Gad c 1 (DAG-WT294)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                         |                                                                        |
|-------------------------|------------------------------------------------------------------------|
| <b>Purity</b>           | > 90%                                                                  |
| <b>Conjugate</b>        | TBD                                                                    |
| <b>Applications</b>     | Immunoassays                                                           |
| <b>Molecular Weight</b> | 13 kDa                                                                 |
| <b>Format</b>           | Liquid                                                                 |
| <b>Concentration</b>    | Batch dependent - please inquire should you have specific requirements |
| <b>Size</b>             | 1 mg                                                                   |
| <b>Buffer</b>           | 20mM Tris, 10mM NaCl, pH8.0                                            |
| <b>Preservative</b>     | None                                                                   |
| <b>Storage</b>          | Store at -80°C. Avoid freeze / thaw cycle                              |

## BACKGROUND

|                 |                                                                            |
|-----------------|----------------------------------------------------------------------------|
| <b>Keywords</b> | <i>Gadus callarias</i> allergen 1; Gad c 1; Beta-parvalbumin; Gad c 1.0101 |
|-----------------|----------------------------------------------------------------------------|