



Recombinant Mangifera indica Man i 1 (DAG-WT295)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Nature	Recombinant
Expression System	E.coli
Purity	> 90%
Conjugate	TBD
Applications	Immunoassays
Molecular Weight	38 kDa
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	20mM Tris-HCl, pH8.0
Preservative	None
Storage	Store at -80°C. Avoid freeze / thaw cycle
Warnings	PLEASE note that this product is intended for research use only; not for diagnostic or clinical use.

BACKGROUND

Keywords	Mangifera indica allergen 1; Man i 1
-----------------	--------------------------------------