



## Bifenthrin [BSA] (DAG004S)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Specificity	N/A
Species	N/A
Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Wet ice

### BACKGROUND

Introduction	Bifenthrin is a pyrethroid insecticide used primarily against the red imported fire ant by influencing its nervous system. It has a high toxicity to aquatic organisms. Bifenthrin is used in agriculture and residential settings, both indoors and outdoors.
Keywords	Bifenthrin;2-Methyl-3-phenylphenyl)methyl (1S,3S)-3-[(Z)-2-chloro-3,3,3-trifluoroprop-1-enyl]-2,2-dimethylcyclopropane-1-carboxylate