



# Levamisole [KLH] (DAG024S)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	N/A
Species	N/A
Conjugate	KLH
Applications	ELISA, LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Wet ice

## BACKGROUND

Introduction	Levamisole ((-)-(S)-2,3,5,6-tetrahydro-6-phenylimidazo (2,1-b) thiazole monohydrochloride) is a veterinary drug used to treat parasitic infestation in animals. It is available as a crystalline white powder, and in pastes, gels, tablets, feed premixes and topical and injectable solutions.
Keywords	Levamisole