



Testosterone [BSA] (DAG3023)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Testosterone, BSA-Conjugated
Species	N/A
Conjugate	BSA
Applications	immunoassay development or other applications.
Format	Liquid
Concentration	>?1.0 mg/mL
Size	1 mg
Buffer	Solution in phosphate buffered saline, pH 7.2 with 0.1% sodium azide.
Preservative	0.1% Sodium Azide
Storage	Store at -20°C.

BACKGROUND

Introduction	Testosterone is a steroid hormone from the androgen group and is found in mammals, reptiles, birds, and other vertebrates. In mammals, testosterone is secreted primarily by the testicles of males and the ovaries of females, although small amounts are also secreted by the adrenal glands. It is the principal male sex hormone and an anabolic steroid.
Keywords	Testosterone