



Recombinant DENV type 3 Envelope Protein [His] (DAG3060)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Dengue Virus Type 3, Envelope Protein (a.a. 281-673), Recombinant. Contains Histidine tag.
Species	DENV
Purity	95% pure (determined by SDS-PAGE). Purified by using immobilized metal-chelate affinity chromatography.
Conjugate	His
Applications	Suitable for use in ELISA, Lateral Flow and Western Blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid
Size	1 mg
Buffer	PBS, pH7.4
Preservative	0.1% Thimerosal, 1 µg/mL of Leupeptin, Aprotinin and Pepstatin A.
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	Dengue virus (DENV) is an enveloped, single-stranded, positive-sense RNA virus that includes four related but distinct serotypes (DENV1, 2, 3, and 4). It encodes three structural proteins (capsid, membrane and envelope) and seven non-structural proteins (NS1, -2a, -2b, -3, -4a, -4b and -5). The envelope (E) glycoprotein mediates virion attachment to the receptor and fusion of
---------------------	--

the virus envelope with the target cell membrane. The recombinant E protein can be used in diagnostic assays for the detection of either primary or secondary dengue infection, to overcome safety issues associated with the use of whole virus.

Keywords	DENV; Dengue Virus Envelope Protein; DENV Envelope Protein; Dengue virus; DENV-3 E Protein; DENV-3
-----------------	--
