



Biotin [HRP] (DAG3197-HRP)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------|---|
| Conjugate | HRP |
| Applications | ELISA |
| Format | Liquid |
| Concentration | Batch dependent - please inquire should you have specific requirements. |
| Size | 1 mg |
| Buffer | 0.01M pH7.4 PBS |
| Preservative | None |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |

BACKGROUND

| | |
|----------|---------------------------|
| Keywords | Biotin;D-biotin;Vitamin H |
|----------|---------------------------|