



## D-Pantothenic Acid [BSA] (DAG3303)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Conjugate	BSA
Applications	LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

Introduction	Pantothenic acid is a member of the class of pantothenic acids that is an amide formed from pantoic acid and beta-alanine. It has a role as a plant metabolite. It is a conjugate acid of a pantothenate.
Keywords	Pantothenic acid; Vitamin B5; D-Pantothenic Acid