



Cloxacillin [HRP] (DAG438S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Cloxacillin, covalently linked to horseradish peroxidase
Nature	Synthetic
Expression System	N/A
Species	N/A
Purity	chromatography, dialysis
Conjugate	HRP
Applications	e.g. allergy test, immunoassays etc.
Molecular Weight	app. 44 kDa
Reconstitution	Reconstitute with 1 ml PBS without azide.
Procedure	None
Format	Lyophilized
Concentration	1 mg/mL working concentration
Size	1 mg
Preservative	None
Storage	Store at -20°C.
Warnings	For research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals.

BACKGROUND

Introduction

Cloxacillin is a semisynthetic antibiotic in the same class as penicillin. Cloxacillin was discovered and developed by Beecham. It is sold under a number of trade names, including Cloxapen, Cloxacap, Tegopen and Orbenin.

Keywords

drug-HRP conjugate, drug-conjugate, Cloxacillin HRP conjugate
