



Cephalosporin [HSA-Biotin] (DAG536S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description Cephalosporin, covalently linked to biotinylated human serum albumin

Species N/A

Purity chromatography, dialysis

Conjugate HSA-Biotin

Applications e.g. allergy test, immunoassays etc.

Molecular Weight app. 65 kDa

Format Liquid

Concentration 1 mg/mL, conjugation ratio is app. 1:10 - 1:40.

Size 1 mg

Buffer PBS

Preservative 0.05% sodium azide

Storage Store at 4°C.

BACKGROUND

Introduction The cephalosporins are a class of β -lactam antibiotics originally derived from Acremonium, which was previously known as "Cephalosporium". Together with cephamycins they constitute a subgroup of β -lactam antibiotics called cephems.

Keywords drug-HSA-Biotin conjugate, drug-conjugate, Cephalosporin HSA-Biotin conjugate