



## Pb<sup>2+</sup> [HRP] (DAGA-005H)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Conjugate	HRP
Applications	ELISA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

Keywords	Lead(2+);Lead (II) ion;Lead(2+) ion;Lead ion;Omaha;Pb <sup>2+</sup>
----------	---