



## Hg<sup>2+</sup> [HRP] (DAGA-007H)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Conjugate	HRP
Applications	ELISA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

Keywords	Mercuric ion;Mercury(2+);Mercury, ion (Hg <sup>2+</sup> )
----------	---