



## Deoxynivalenol [OVA] (DAGA-0120)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Conjugate</b>	OVA
<b>Applications</b>	ELISA
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.