



Sterigmatocystin [HRP] (DAGA-015H)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|-------------------|---|
| Nature | Synthetic |
| Expression System | N/A |
| Conjugate | HRP |
| Applications | ELISA |
| Procedure | None |
| Format | Liquid |
| Concentration | Batch dependent - please inquire should you have specific requirements. |
| Size | 1 mg |
| Buffer | 0.01M pH7.4 PBS |
| Preservative | None |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |

BACKGROUND

| | |
|----------|------------------|
| Keywords | Sterigmatocystin |
|----------|------------------|