



Recombinant DENV type 1 polyprotein protein (DENV1/CN/GZ13/2014, Fragment) [His] (DAGA-1043)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Recombinant 6xHis tagged DENV type 1 polyprotein protein (DENV1/CN/GZ13/2014, Fragment) (UniProtKB No. A0A0C5CHG7) was expressed.
Nature	Recombinant
Expression System	E. coli
Purity	>85%
Conjugate	His
Applications	ELISA, Immunogen
Recommended Usage	This item requires custom production and lead time is between 4-8 weeks. We can custom produce based on your specifications.
Procedure	None
Format	Lyophilized
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg, 5 mg
Buffer	Each vial contains 100 µg of purified protein in PBS
Preservative	None
Storage	Store at 2-8 °C short term and -20°C long term. It is recommended to aliquot needed volumes

and store at -20°C. Avoid repeated freeze/thaw cycles.

BACKGROUND

Keywords DENV type 1;Dengue Virus type 1;Dengue Virus 1+2+3+4;Dengue Virus Type 1, 2, 3, 4;Group IV;Flaviviridae;Flavivirus;Dengue virus;Dengue Virus;DENV

ANTIGEN GENE INFORMATION

UniProt ID A0A0C5CHG7
