



# Recombinant DENV type 3 Nonstructural protein 5 (Fragment) [His] (DAGA-1059)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	Recombinant 6xHis tagged DENV type 3 Nonstructural protein 5 (Fragment) (UniProtKB No. A0A0K0NLD9) was expressed.
<b>Purity</b>	>85%
<b>Conjugate</b>	His
<b>Applications</b>	ELISA, Immunogen
<b>Recommended Usage</b>	This item requires custom production and lead time is between 4-8 weeks. We can custom produce based on your specifications.
<b>Format</b>	Lyophilized
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg, 5 mg
<b>Buffer</b>	Each vial contains 100 µg of purified protein in PBS
<b>Preservative</b>	None
<b>Storage</b>	Store at 2-8 °C short term and -20°C long term. It is recommended to aliquot needed volumes and store at -20°C. Avoid repeated freeze/thaw cycles.

## BACKGROUND

<b>Keywords</b>	DENV type 3; Dengue Virus type 3; Dengue Virus 1+2+3+4; Dengue Virus Type 1, 2, 3, 4; Group IV; Flaviviridae; Flavivirus; Dengue virus; Dengue Virus; DENV
-----------------	--

# GENE INFORMATION

UniProt ID                      A0A0K0NLD9

---