



## Dimethyl Phthalate [KLH] (DAGA-108K)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|               |   |
|---------------|---|
| Conjugate     | KLH   |
| Applications  | Immunogen   |
| Format        | Liquid  |
| Concentration | Batch dependent - please inquire should you have specific requirements.               |
| Size          | 1 mg  |
| Buffer        | 0.01M pH7.4 PBS   |
| Preservative  | None  |
| Storage       | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |

### BACKGROUND

|          |                    |
|----------|--------------------|
| Keywords | Dimethyl phthalate |
|----------|--------------------|