



Vitamin B1 [KLH] (DAGA-120K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Also known as thiamine, vitamin B1 plays an important role in helping the body convert carbohydrates and fat into energy.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Aberic acid; Thiamin; Thiamine; Vitamin B1;
----------	---