



Vitamin B3 [KLH] (DAGA-122K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Niacin, also known as vitamin B3 and nicotinic acid, is an organic compound with the formula C ₆ H ₅ NO ₂ and, depending on the definition used, one of the 20 to 80 essential human nutrients. Pharmaceutical and supplemental niacin are primarily used to treat hypercholesterolemia (high cholesterol) and pellagra (niacin deficiency).
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords Niacin; Nicotinamide; Nicotinic acid;