



Vitamin B5 [KLH] (DAGA-123K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Pantothenic acid; pantothenate; vitamin B5; Pantothenic; Ca pantothenaat; kyselinapantothenova; d-Pantothenic Acid; Vitamin B5
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	D-pantothenic acid; Pantothenic acid; Vitamin B5;
----------	---