



## Vitamin C [KLH] (DAGA-125K)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Antigen Description	Vitamin C (L ascorbic acid and its reduced form, dehydroascorbic acid) is a water-soluble nutrient and vitamin essential for life and for maintaining optimal health.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

**Keywords** 3 oxo L gulofuranolactone; Ascorbic acid; Dehydroascorbic acid; Salt ascorbate; Vitamin C;